

Sarvin Sanat Mana Company

Supplier Sealing for Oil and Gas / Melting Metals



Ceramic 99% / SGL

The Corosion Test For Ceramic in Reagent		
Test Enviroment		Ceramic 99%
Reagent Concentration	Temperature (C)	Ceramic
98% H2SO4	100	65
50% NaOH	100	75
53% HF	25	20
85% H3PO4	100	>1000
70% HNO3	100	7
45% KOH	100	60
25% HCL	70	72
10% HF + 57% HNO3	25	16